



Sizing Chart

Body Systems

TheraTogs are designed to fit typical body types within a given weight range. If either length measurement exceeds the maximum shown for wearer's weight range, select the next larger size. For example, for a child weighing 40 pounds with a Hip-to-Knee measurement of 13", select the #300 Pediatric model.

Series Size No. & Name		U.S. Units			Metric Units		
		Weight (Lbs)	Max Shoulder-to-Hip (in)	Max Hip-to-Knee Joint (in)	Weight (Kilos)	Max Shoulder-to-Hip (cm)	Max Hip-to-Knee Joint (cm)
50	Small Infant	0 - 15	8½	7½	0 - 6.8	21.6	19.1
100	Infant	15 - 25	10½	9¾	6.80 - 11.3	26.7	24.8
200	Pre-School	25 - 45	13½	11¾	11.3 - 20.4	34.3	29.9
300	Pediatric	45 - 65	15	18½	20.4 - 29.4	38.2	47.0
400	Petite Adult	65 - 110	17½	20¾	29.4 - 49.9	44.5	52.8
600	Adult Female Small	110 - 140	20½	23¾	49.9 - 63.5	52.1	60.4
610	Adult Female Med	140 - 170	20½	23¾	63.5 - 77.1	52.1	60.4
620	Adult Female Large	170 - 200	22¾	26¼	77.1 - 90.9	57.8	66.7
700	Adult Male Small	110 - 145	19¼	23¾	49.9 - 65.8	48.9	60.4
710	Adult Male Med	145 - 170	19¼	23¾	65.8 - 77.1	48.9	60.4
720	Adult Male Large	170 - 190	23¼	26¼	77.1 - 86.1	59.1	66.7
730	Adult Male X Large	190 - 240	22½	26¼	86.1 - 108.9	57.2	66.7

The Hip measurement point is the ASIS -- the "front hipbone" (anterior superior iliac spine – the front top of the iliac crest). The Knee measurement point is the top of the patella ("knee cap").

Gait Trainer Shoe Size Chart

Use these shoe size ranges to determine the best model selection for a Gait Trainer System:

Model #	Description	US Size	Europe	UK	Inches	CM
#204	Toddler	3 - 8.5	18 - 25	2 - 7.5	4.125 - 5.75	10.5 - 14.6
#304	Toddler	9 - 13.5	25 - 31	8 - 12.5	6.125 - 7.625	15.6 - 19.4
	Girls & Boys	1 - 2	32 - 33	13 - 1	7.75 - 8.125	19.7 - 20.6
#314	Girls & Boys	3 - 6	34 - 38	2 - 5	8.5 - 9.5	21.6 - 24.1
#404	Girls/Women	5 - 8.5	35 - 39		8.5 - 9.6875	
	Boys/Men	3 - 6.5	35 - 39		8.5 - 9.5	
#614	Women	9 - 12	42 - 45		10.25 - 11.25	
	Men	7 - 9.5	40 - 43		9.6 - 10.4	
#624	Men	10 - 15	43 - 48		10.5 - 12.2	

THERATOGS™ IS A TRADEMARK OF THERATOGS, INC. ALL RIGHTS RESERVED. U.S. AND FOREIGN PATENTS PENDING.

TheraTogs, Inc • 305 Society Drive, Unit C-3 • Telluride, CO 81435
 (888) 634-0495 • (970) 728-7078 • (970) 728-7028 fax • www.theratogs.com • info@theratogs.com