




TheraTogs Technical Highlight

August 2007

Volume 1, Issue 1

"I have had them on my son for 6 weeks now. Bottom line is I LOVE the new straps...they look and feel and perform like new!"

- Jill, caregiver



TogRite Strapping features a unique honeycomb pattern to maximize grip and breathability.

www.theratogs.com

PHONE:
1-888-634-0495

FAX:
(970) 728-7028

E-MAIL:
info@theratogs.com

Introducing TheraTogs *TogRite* Strapping

TheraTogs, Inc. is pleased to announce the introduction of **TogRite Strapping** to the TheraTogs product line.

The increasing demand and expanding range of rehabilitation applications for TheraTogs systems have created a need for a more robust strapping system.

This new and improved strapping will replace the

current MediWrap material (the stretchy white strapping) used in most TheraTogs systems.

TogRite strapping is comprised of much stronger materials, which expands the ability for TheraTogs applications to accommodate larger client sizes, and to support limb and joint positioning that requires more force.

TogRite strapping also has increased tolerance to a wider range of

conditions, including machine washing.

This technical newsletter explains the properties and advantages of this new upgrade.

The introduction of the new strapping will initially be included in Pre-school through Adult sized systems for all two-inch (2") and three-inch (3") wide strapping.

TogRite Technical Performance

TogRite is a latex free "external muscle" material that incorporates a combination of nylon and spandex with a backing of inert silicone grip surface. This combination of components is designed to resemble a sarcomere structure and the gripping properties of the silicone provide a unique replica of a functioning muscle.

Standardized ASTM tests comparing the new TogRite product to the current MediWrap strapping – as well as other strapping materials on the market – demonstrate that when placed under the same load, the TogRite product maintains longitudinal position and rebounds to the original length. Even after extended use, TogRite strapping shows no "stretch out" or loss of strength and resistance.

This is a significant improvement from previous strapping materials – obviously the stronger and more resilient the strapping, the more effective the TheraTogs application will be for its wearer.

Additional field tests demonstrated the benefits of the special weave of the external surface of the TogRite materials. Performance of the strap remained at 100% rebound, even when the material was purposefully frayed or cut.

This unique combination of materials provides assurance to the therapist that the strap will remain in the correct placement.

These straps are so robust that strap length requirements can be decreased by up to 20% when compared to previous strap sizes provide in TheraTogs systems.



Adult and larger limb rotational requirements have changed the technical requirements for our strapping. TogRite is fantastic for adult patients with gait issues.



TogRite Product Characteristics

- Machine Washable
- Reduced Fraying and Fuzzing
- Packaged in Rolls to allow YOU to Customize Lengths
- Longer Lasting -- new material construction!
- More Consistent Stretch & Rebound
- Stronger Resistance
- Improved Grip via New Inert Silicone Backing

TheraTogs TogRite Replacement Policy

White strapping in TheraTogs systems purchased between July, 2006 and December, 2007 exhibiting over stretching, lack of rebound or in-sufficient rotational tension is available for complimentary replacement. Please refer to www.theratogs.com/warranty, Warranty/Return Policy for additional details.

TogRite Availability in Rolls

The new TogRite rolls will be available to customers starting August 2007.

All new TheraTogs systems purchased (Pres-School through Adult) will include the new two inch and three inch wide TogRite strapping.

To allow for therapist flexibility during treatment, the TogRite straps provided in systems will

now be packaged in rolls rather than cut to specific length.

TogRite will also be available (in rolls) for purchase by clinicians to replenish strapping in a Clinical Sizing Kit, to replace client strapping that has lost effectiveness, or to address TheraTogs applications needing more resilient strapping materials.



Introductory pricing for these rolls is available through the balance of 2007.

2" Wide x 1 Yard \$25

3" Wide x 1 Yard \$36

2" Wide x 3 Yards \$62

3" Wide x 3 Yards \$95